



CareElements

Showering

Nursing bench

In situations where a person needs to be bathed while lying down, it is important that the nursing bench can be lowered to a suitable height for transfer, for example from a wheelchair, and then raised to a suitable working height for the carer. It is equally important to ensure that the bench is sufficiently comfortable and safe.

The electrically powered nursing bench can be adjusted by 700 mm simply by using the remote control.

Comfortable surfaces – the bench, the end caps and the safety rail brackets are covered with polyurethane foam.

The slightly curved form makes the bench extremely comfortable and safe.

The water collection tray is equipped with an integrated water outlet that prevents water from splashing on carers or the floor.

The integrated gas cylinder enables the bench to be folded easily against the wall when not in use to maximise floor space.




Winner
The Danish Design Prize 2007


① **Manually adjustable back support** at both ends for increased comfort.

② **The flexible hose** attached to the centre mounted outlet of the water collection tray can be connected to either an outlet in the floor or on the wall.


③ **The safety rail** can be folded easily up or down.



Mounting: Electric model
Fixing bolts for solid walls are included.



Mounting: Fixed model
Fixing bolts for solid walls are included.



Materials: Aluminium, polyurethane foam with marine ply core, powder lacquered steel

Maximum weight: 200 kg

Instructions for mounting products onto different wall types can be downloaded from www.pressalitcare.com

Nursing bench – adjustable height


- 1800 mm, 700 mm height adj., electric motor **R8518**
- 1300 mm, 700 mm height adj., electric motor **R8513**
- 1800 mm, 700 mm height adj., el. motor, safety rail **R8538**
- 1300 mm, 700 mm height adj., el. motor, safety rail **R8533**

When mounted at the recommended height, the nursing bench's height can be adjusted from 300 to 1000 mm above floor.

If safety rail is to be pre-mounted, please be sure to indicate the correct item number when ordering.

Incl. remote control and flexible hose, 1500 mm.

For placement of electrical installations, see page 68.



R8518

Projection when folded:
With safety rail 330 mm
Without safety rail 300 mm



R8518


Nursing bench – fixed height

- 1800 mm **R8508**
- 1300 mm **R8503**
- 1800 mm, safety rail **R8528**
- 1300 mm, safety rail **R8523**

When the bench is mounted at the recommended height, the lying surface height is 865 mm above floor.

If safety rail is to be pre-mounted, please be sure to indicate the correct item number when ordering.

Incl. flexible hose, 1500 mm.



R8523


Projection when folded:
With safety rail 310 mm
Without safety rail 280 mm




R8508

Safety rail


- 1636 mm **R8511**
- 1200 mm **R8510**



R8511




R8510



R8511


Accessories for nursing bench

RT610




Shower head:
See page 40.

RT611




Shower mixer:
See page 40.

RT640




RT641




Rack system:
Racks, shelves, boxes and knobs. See page 48

RT9xx



RT7xx



Baskets:
See page 50.